

VITANIX CENTRIFUGE VC-60



Benchtop Automatic Centrifuge

VC-60

OVERVIEW

The centrifuge is perfect for medical, hospital, pathology and institutional laboratories together with a broad well-chosen of rotors with adaptors.

Features & Benefits

- Microprocessor controlled with built in speedometer & speed control, Safety Interlock
- Digital Display & control of Timing Range 0~99min
- Brushless motor.
- Automatic balance. & rotor Imbalance Detetction
- Low noise and low temperature-increase in chamber.
- Metal external case.
- Maintenance free induction drive motor.
- Rubber suspensions reduce the vibrations.
- Low speed centrifuge with wide choice of rotors.
- Extremely quiet operation.
- Door lid motorized.
- Automatic door opening through gas spring.

Specification

MAX.SPEED	6000 rpm	DISPLAY	LCD
MOTOR	BRUSHLESS MOTOR	POWER	AC220V, 50/60Hz,5A
TIME RANGE	0~99min	WEIGHT	20kg
NOISE	<55db	DIMENSION	450x365x280mm(LxWxH)
MAX RCF	5080xg	Acceleration/ Deceleration Rates	1--10

R o t o r L i s t

ROTOR	DESCRIPTION	MAX.SPEED	MAX RCF
ANGLE ROTOR	12X10ml	5000rpm	2095xg
	18X10ml	5000rpm	2240xg
	12X15ml	6000rpm	5080xg
	12X20ml	5000rpm	2200xg
	6X50ml	5000rpm	4745xg
	4X50ml	5000rpm	4745xg
	24X15/10/5ml	4000 /min	2160 xg
SWING BUCKET ROTOR	10X10ml	5000rpm	2500xg
	12X10ml	5000rpm	2600xg
	16X10ml	5000rpm	2600xg
	8X15ml	5000rpm	2690xg
	4X50ml	5000rpm	2420xg



VITANIX

future
of
healthcare



+38664150661 | INFO@VITANIX.SI | WWW.VITANIX.SI

ADDRESS: Mlinska ulica 1 2000 Maribor Slovenia Europe

VITANIX D.o.o reserves the right to make any changes to the information, including specifications, without notice

